

Product Datasheet

SUPT16H antibody (orb19521)

Description

Goat polyclonal antibody to SUPT16H



Species/Host

Goat

Reactivity

Bovine, Human, Mouse, Porcine

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

orb19521 (0.1ug/ml)
staining of HeLa cel...

Target

FACT / SUPT16H

Preservatives

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.

Storage


Aliquot and store at -20°C. Minimize freezing and thawing.

Note

For research use only

Application notes

ELISA: Peptide ELISA: antibody detection limit dilution 1:8000. WB: Approx 150kDa band observed in lysates of cell line HeLa (calculated MW of 120kDa according to NP_009123.1). Recommended concentration: 0.1-0.3 µg/ml. An additional band of unknown identity was also consistently observed at 30kDa. This band was successfully blocked by incubation with the immunising peptide. We would appreciate any feedback from people in the field - have any such results been reported with other antibodies/lysates Have any further splice variants/modified forms been reported



Immunohistochemical
staining of Human
Ki...

Protein Sequence

DLTNKEGKKPEEK

Clonality

Polyclonal

MW

120

Entrez

11198

Dilution Range

ELISA: 1:8000, WB: 0.1-0.3 µg/ml

Expiration Date

12 months from date of receipt.