

## Product Datasheet

# Human Carbonic anhydrase-related protein 11 (CA11) ELISA Kit (orb1952290)

|                            |  |
|----------------------------|--|
| <b>Description</b>         | Human Carbonic anhydrase-related protein 11 (CA11) ELISA Kit   |
| <b>Tested Applications</b> | ELISA  |
| <b>Note</b>                | For research use only  |
| <b>Application notes</b>   | Detect Range: 0.78-50 U/mL. Sensitivity: 0.29 U/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 50 µL. Sample Type: serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids. |
| <b>Expiration Date</b>     | Please enquire.  |