

## Product Datasheet

# Human Interferon Gamma (IFN $\gamma$ ) ELISA Kit (orb1952292)

<b>Description</b>	Human Interferon Gamma (IFN $\gamma$ ) ELISA Kit
<b>Tested Applications</b>	ELISA
<b>Note</b>	For research use only
<b>Application notes</b>	Detect Range: 31.2-2000 pg/mL. Sensitivity: 11.74 pg/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 50 $\mu$ L. Sample Type: serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
<b>Expiration Date</b>	Please enquire.