

Product Datasheet

Rat Milk Fat Globule EGF Factor 8 (MFGE8) ELISA Kit (orb1952305)

Description	Rat Milk Fat Globule EGF Factor 8 (MFGE8) ELISA Kit
Reactivity	Rat
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 31.2-2000 pg/mL. Sensitivity: 12.8 pg/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 50 µL. Sample Type: serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Expiration Date	Please enquire.