

## Product Datasheet

### Rat Decorin (DCN) ELISA Kit (orb1952332)

<b>Description</b>	Rat Decorin (DCN) ELISA Kit
<b>Reactivity</b>	Rat
<b>Tested Applications</b>	ELISA
<b>Note</b>	For research use only
<b>Application notes</b>	Detect Range: 0.156-10 ng/mL. Sensitivity: 0.62 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 50 µL. Sample Type: serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
<b>Expiration Date</b>	Please enquire.