

Product Datasheet

Mouse NACHT, LRR and PYD Domains-Containing Protein 3(NLRP3) ELISA Kit (orb1952361)

Description	Mouse NACHT, LRR and PYD Domains-Containing Protein 3(NLRP3) ELISA Kit
Reactivity	Mouse
Tested Applications	ELISA
Note	For research use only
Application notes	Detect Range: 0.156-10 ng/mL. Sensitivity: 0.055 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 50 µL. Sample Type: serum,plasma,tissuehomogenates,celllysates,cellculturesupernatesorotherbiologicalfluids.
Expiration Date	Please enquire.