

Product Datasheet

Human Soluble terminal complement complex C5b-9 (C5b-9) ELISA Kit (orb1952393)

| | |
|----------------------------|--|
| Description | Human Soluble terminal complement complex C5b-9 (C5b-9) ELISA Kit |
| Reactivity | Human |
| Tested Applications | ELISA |
| Note | For research use only |
| Application notes | Detect Range: 31.2-2000 ng/mL. Sensitivity: 12.6 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 50 µL. Sample Type: serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids. |
| Expiration Date | Please enquire. |