

## Product Datasheet

### Mouse Defensin-6(DEFA6) ELISA kit (orb1952478)

<b>Description</b>	Mouse Defensin-6(DEFA6) ELISA kit
<b>Reactivity</b>	Mouse
<b>Tested Applications</b>	ELISA
<b>Note</b>	For research use only
<b>Application notes</b>	Detect Range: 0.312-20 ng/mL. Sensitivity: 0.126 ng/mL. Assay Time: 1-4.5 h. Detection wavelength: 450 nm. Sample Volume: 50 µL. Sample Type: serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
<b>Expiration Date</b>	Please enquire.