

Product Datasheet

Recombinant Human Otoraplin (orb1952538)

Description	Recombinant Human Otoraplin
Endotoxins	Less than 1 EU/ug of rHuOTOR as determined by LAL method.
Form/Appearance	Lyophilized from a 0.2 µm filtered concentrated solution in 20 mM PB, pH 7.4, 150mM NaCl.
Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles.- 12 months from date of receipt, -20 to -70 °C as supplied.- 1 month, 2 to 8 °C under sterile conditions after reconstitution.- 3 months, -20 to -70 °C under sterile conditions after reconstitution.
Note	For research use only
Protein Sequence	VHGIFMDRLA SKKLCADDEC VYTISLASAQ EDYNAPDCRF INVKKGQIQY VYSKLVKENG AGEFWAGSVY GDGQDEMGVV GYFPRNLVKE QRVYQEATKE VPTTDIDFFC E
Purity	> 97 % by SDS-PAGE and HPLC analyses.
MW	Approximately 12.7 kDa, a single non-glycosylated polypeptide chain containing 111 amino acids.
Expiration Date	6 months from date of receipt.