

## Product Datasheet

# Recombinant Bovine Enterokinase Light Chain, His (orb1952554)

<b>Description</b>	Recombinant Bovine Enterokinase Light Chain, His
<b>Endotoxins</b>	Less than 1 EU/ug of rBoEKL, His as determined by LAL method.
<b>Form/Appearance</b>	50 mM Tris-HCl, pH 8.0, 0.5 M NaCl and 50 % glycerol.
<b>Storage</b>	Use a manual defrost freezer and avoid repeated freeze-thaw cycles.- 6 months from date of receipt, -20 to -70 °C as supplied.- 3 months, -20 to -70 °C under sterile conditions after opening.
<b>Note</b>	For research use only
<b>MW</b>	Approximately 28 kDa, a single non-glycosylated polypeptide chain containing 241 amino acids with 6 x His at C-terminus.
<b>Expiration Date</b>	6 months from date of receipt.