

Product Datasheet

Recombinant Bovine Enterokinase Light Chain, Yeast (orb1952555)

Description	Recombinant Bovine Enterokinase Light Chain, Yeast
Endotoxins	Less than 1 EU/ug of rBoEKL, His as determined by LAL method.
Form/Appearance	50 mM Tris-HCl, pH 8.0, 0.5 M NaCl and 50 % glycerol.
Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles.- 6 months from date of receipt, -20 to -70 °C as supplied.- 3 months, -20 to -70 °C under sterile conditions after opening.
Note	For research use only
MW	Approximately 43.0 kDa, a single glycosylated polypeptide chain containing 235 amino acids.
Expiration Date	6 months from date of receipt.