
Product Datasheet

Human Klotho beta / KLB Protein, His Tag (orb1952613)

Description

Human Klotho beta Protein, His Tag is

(HEK293). It contains AA Met 30...

Reactivity

Human

Endotoxins

1.0 EU per µg

Conjugation

Unconjugated

Target

Klotho beta

Preservatives

Powder

Form/Appearance

Powder

Storage

-20°C to -70°C for 12 months in lyophilized state;. -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower

Tag

C-10×His

Note

For research use only

Application notes

This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 112.9 kDa. The protein migrates as 65 kDa,95 kDa and 130 kDa when calibrated against Star Ribbon Pre-stained Protein Marker under reducing (R) condition (SDS-PAGE).

Purity

90%

Source

Human Klotho beta Protein, His Tag is expressed from human 293 cells (HEK293). It contains AA Met 30 - Leu 997 (Accession # Q86Z14-1).

MW

112.9 kDa

Uniprot ID

Q86Z14

Expiration Date

6 months from date of receipt.