
Product Datasheet

**Cynomolgus PVRIG Protein, His Tag (MALS
verified) (orb1952618)**

Description

Cynomolgus PVRIG Protein, His Tag

(HEK293). It contains AA Thr 41 ...

Reactivity

Monkey

Endotoxins

0.1 EU per µg

Conjugation

Unconjugated

Target

PVRIG

Preservatives

Powder

Form/Appearance

Powder

Storage

-20°C to -70°C for 12 months in lyophilized state;. -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower

Tag

C-10×His

Note

For research use only

Application notes

This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 15.6 kDa. The protein migrates as 20-27 kDa when calibrated against Star Ribbon Pre-stained Protein Marker under reducing (R) condition (SDS-PAGE) due to glycosylation.

Purity

90%

Source

Cynomolgus PVRIG Protein, His Tag is expressed from human 293 cells (HEK293). It contains AA Thr 41 - Asp 171 (Accession # A0A2K5WVV8-1).

MW

15.6 kDa

Uniprot ID

[A0A2K5WVV8](#)

Expiration Date

6 months from date of receipt.