

---

## Product Datasheet

**Rhesus macaque CD155 / PVR Protein, Fc Tag  
(MALS verified) (orb1952619)**

**Description**

Rhesus macaque CD155 Protein,

293 cells (HEK293). It contains AA Asp ...

**Reactivity**

Monkey

**Endotoxins**

1.0 EU per µg

**Conjugation**

Unconjugated

**Target**

CD155

**Preservatives**

Powder

**Form/Appearance**

Powder

**Storage**

-20°C to -70°C for 12 months in lyophilized state;. -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower

**Tag**

C-Fc

**Note**

For research use only

**Application notes**

This protein carries a human IgG1 Fc tag at the C-terminus. The protein has a calculated MW of 60.8 kDa. The protein migrates as 70-90 kDa when calibrated against Star Ribbon Pre-stained Protein Marker under reducing (R) condition (SDS-PAGE) due to glycosylation.

**Purity**

90%

**Source**

Rhesus macaque CD155 Protein, Fc Tag is expressed from human 293 cells (HEK293). It contains AA Asp 28 - Asn 343 (Accession # Q0MSE6-1).

**MW**

60.8 kDa

**Uniprot ID**
**Q0MSE6**
**Expiration Date**

6 months from date of receipt.