

---

## Product Datasheet

**Human Kallikrein 3 Protein, His Tag (active enzyme) (orb1952623)**

**Description**

Human Kallikrein 3 Protein, His Tag

(HEK293).

**Reactivity**

Human

**Endotoxins**

1.0 EU per µg

**Conjugation**

Unconjugated

**Target**

Kallikrein 3

**Preservatives**

Powder

**Form/Appearance**

Powder

**Storage**

-20°C to -70°C for 12 months in lyophilized state;. -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower

**Tag**

C-10×His

**Note**

For research use only

**Application notes**

This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 28.7 kDa. The protein migrates as 33-36 kDa when calibrated against Star Ribbon Pre-stained Protein Marker under reducing (R) condition (SDS-PAGE) due to glycosylation.

**Purity**

95%

**Source**

Human Kallikrein 3 Protein, His Tag is expressed from human 293 cells (HEK293).

**MW**

28.7 kDa

**Uniprot ID****P07288****Expiration Date**

6 months from date of receipt.