
Product Datasheet

Human IL-3 Protein, premium grade (orb1952634)

Description

Human IL-3 Protein, premium grade is

AA Ala 20 - Phe 152 (...

Reactivity

Human

Endotoxins

0.01 EU per µg

Conjugation

Unconjugated

Target

IL-3

Preservatives

Powder

Form/Appearance

Powder

Storage

-20°C to -70°C for 24 months in lyophilized state;. -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower

Tag

Native

Note

For research use only

Application notes

This protein carries no "tag". The protein has a calculated MW of 15.1 kDa. The protein migrates as 14 kDa±3 kDa when calibrated against Star Ribbon Pre-stained Protein Marker under reducing (R) condition (SDS-PAGE).

Purity

95%

Source

Human IL-3 Protein, premium grade is expressed from E. coli cells. It contains AA Ala 20 - Phe 152 (Accession # P08700-1).

MW

15.1 kDa

Uniprot ID
P08700
Expiration Date

6 months from date of receipt.