

Product Datasheet

GMP Human Laminin 521 Protein (orb1952636)

Description	GMP Human Laminin 521 Protein is expressed from human 293 cells (HEK293).
Reactivity	Human
Endotoxins	10 EU/mg
Conjugation	Unconjugated
Target	Laminin 521
Preservatives	Powder-GMP
Form/Appearance	Powder-GMP
Storage	-20°C to -70°C for 5 years in lyophilized state;. -70°C for 12 months under sterile conditions after reconstitution. Upon receipt, store it immediately at -20°C or lower for long term storage
Note	For research use only
Application notes	This protein .
Purity	95%
Source	GMP Human Laminin 521 Protein is expressed from human 293 cells (HEK293).
Uniprot ID	P55268, P11047, O15230
Expiration Date	6 months from date of receipt.