

## Product Datasheet

### Recombinant Human RIP-3 (orb1952676)

<b>Description</b>	Recombinant Human RIP-3
<b>Reactivity</b>	Human
<b>Tested Applications</b>	ELISA, WB
<b>Preservatives</b>	Lyophilized from a 0.2 µm filtered solution in 10 mM Hepes, 500 mM NaCl with 5% trehalose, pH 7.4
<b>Storage</b>	The lyophilized protein is stable at -20°C for up to 1 year, the protein solution can be stored at 2-8°C for up to 1 week. For extended storage, it is recommended to further dilute in a buffer containing 0.1% BSA (make sure BSA will NOT interfere with your experiment ) and store in working aliquots at -20°C to -80°C for 6 months. Avoid repeated freeze/thaw cycle.
<b>Tag</b>	N-terminal His Tag
<b>Note</b>	For research use only
<b>Application notes</b>	Centrifuge the vial prior to opening, reconstitute in sterile distilled water or buffer of interest to concentration of 0.1-1mg/ml by gently pipetting 2-3 times, don't vortex
<b>Purity</b>	Greater than 95% as determined by SDS-PAGE
<b>Source</b>	E.Coli
<b>MW</b>	33.2 kDa
<b>Uniprot ID</b>	<a href="https://www.uniprot.org/entry/Q9Y572">Q9Y572</a>
<b>Expiration Date</b>	6 months from date of receipt.