

Product Datasheet

Rabbit anti-PEG (Polyethylene glycol) IgM ELISA Kit (orb1952751)

Description	Rabbit anti-PEG (Polyethylene glycol) IgM ELISA Kit
Reactivity	Rabbit
Range	0.156-10ng/ml
Target	anti-PEG IgM
Storage	2-8 °C for 6 months
Note	For research use only
Sample Types	Serum, Plasma, Cell Culture Supernatant, cell or tissue lysate, Other liquid samples
Sensitivity	0.094ng/ml
Expiration Date	Please enquire.

