

Product Datasheet

Porcine HS-INS(High sensitive Insulin)Accquant ELISA Kit (orb1952781)

Description	Porcine HS-INS(High sensitive Insulin)Accquant ELISA Kit
Reactivity	Rat
Range	7.813-500pg/ml
Target	INS
Storage	2-8 °C for 6 months
Note	For research use only
Sample Types	Serum, Plasma, Cell Culture Supernatant, cell or tissue lysate, Other liquid samples
Uniprot ID	P01322
Sensitivity	4.688pg/ml
Expiration Date	Please enquire.

