

Product Datasheet

Rat OC/BGP (Osteocalcin) ELISA Kit (orb1952784)

Description	Rat OC/BGP (Osteocalcin) ELISA Kit
Reactivity	Rat
Range	0.78-50pg/ml
Target	OC/BGP
Storage	2-8 °C for 6 months
Note	For research use only
Sample Types	Serum, Plasma, Cell Culture Supernatant, cell or tissue lysate, Other liquid samples
Uniprot ID	P04640
Sensitivity	0.469pg/ml
Expiration Date	Please enquire.