

Product Datasheet

Rat MCP-1 (Monocyte Chemotactic Protein 1) ELISA Kit (orb1952791)

Description	Rat MCP-1 (Monocyte Chemotactic Protein 1) ELISA Kit
Reactivity	Rat
Range	3.125-200pg/ml
Target	MCP-1
Storage	2-8 °C for 6 months
Note	For research use only
Sample Types	Serum, Plasma, Cell Culture Supernatant, cell or tissue lysate, Other liquid samples
Uniprot ID	P14844
Sensitivity	1.875pg/ml
Expiration Date	Please enquire.