

Product Datasheet

Mouse Ldhb (L-lactate dehydrogenase B chain) ELISA Kit (orb1952842)

Description	Mouse Ldhb (L-lactate dehydrogenase B chain) ELISA Kit
Reactivity	Mouse
Range	0.781-50ng/ml
Target	Ldhb
Storage	2-8 °C for 6 months
Note	For research use only
Sample Types	Serum, Plasma, Cell Culture Supernatant, cell or tissue lysate, Other liquid samples
Uniprot ID	P16125
Sensitivity	0.469ng/ml
Expiration Date	Please enquire.