

## Product Datasheet

# Human VEGFR1 (Vascular Endothelial Growth Factor Receptor 1) ELISA Kit (orb1952940)

<b>Description</b>	Human VEGFR1 (Vascular Endothelial Growth Factor Receptor 1) ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	31.25-2000pg/ml
<b>Target</b>	VEGFR1
<b>Storage</b>	2-8 °C for 6 months
<b>Note</b>	For research use only
<b>Uniprot ID</b>	<a href="#">P17948</a>
<b>Sensitivity</b>	18.75pg/ml
<b>Expiration Date</b>	Please enquire.