

## Product Datasheet

# Human IGF-2 (Insulin-like growth factor II) ELISA Kit (orb1952942)

<b>Description</b>	Human IGF-2 (Insulin-like growth factor II) ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	31.25-2000pg/ml
<b>Target</b>	IGF-2
<b>Storage</b>	2-8 °C for 6 months
<b>Note</b>	For research use only
<b>Uniprot ID</b>	<b>P01344</b>
<b>Sensitivity</b>	18.75pg/ml
<b>Expiration Date</b>	Please enquire.