

Product Datasheet

Human TNFRSF9/4-1BB (Tumor necrosis factor receptor superfamily member 9) ELISA Kit (orb1952999)

Description	Human TNFRSF9/4-1BB (Tumor necrosis factor receptor superfamily member 9) ELISA Kit
Reactivity	Human
Range	9.375-600pg/ml
Target	TNFRSF9
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	Q07011
Sensitivity	5.625pg/ml
Expiration Date	Please enquire.