

Product Datasheet

IgG1 Fc Protein, Mouse, Recombinant (C102S) (orb1953346)

Description	IgG1 Fc Protein, Mouse, Recombinant (C102S) is expressed in HEK293 mammalian cells. The predicted molecular weight is 25.8 kDa and the accession number is AAK53870.1.
Storage	-20°C
Tag	Tag Free
Note	For research use only
Application notes	A Certificate of Analysis (CoA) containing reconstitution instructions is included with the products. Please refer to the CoA for detailed information.
Purity	98.00%
MW	25.8 kDa (predicted); 32 kDa (reducing condition, due to glycosylation)
Uniprot ID	AAK53870.1
Expiration Date	6 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)