

Product Datasheet

Monkeypox Virus (MPXV) Envelope protein E8L Protein (His) (orb1953387)

Description	Monkeypox Virus (MPXV) Envelope protein E8L Protein (His) is expressed in Baculovirus insect cells with His tag. The predicted molecular weight is 35.75 kDa and the accession number is QNP13700.1.
Storage	-20°C
Tag	C-His
Note	For research use only
Application notes	A Certificate of Analysis (CoA) containing reconstitution instructions is included with the products. Please refer to the CoA for detailed information.
Purity	98.00%
MW	35.75 kDa (predicted); 67.64 kDa and 37.35 kDa (reducing conditions)
Uniprot ID	P13700
Expiration Date	6 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)