

Product Datasheet

Cathepsin H Protein, Human, Recombinant (His) (orb1953464)

Description	Proinflammatory cytokines decreased the expression of CTSH in human islets and primary rat beta-cells, and overexpression of CTSH protected insulin-secreting cells against cytokine-induced apoptosis. Cathepsin H Protein, Human, Recombinant (His) is expressed in HEK293 mammalian cells with His tag. The predicted molecular weight is 36.62 kDa and the accession number is NP_004381.2.
Storage	-20°C
Tag	C-His
Note	For research use only
Application notes	A Certificate of Analysis (CoA) containing reconstitution instructions is included with the products. Please refer to the CoA for detailed information.
Purity	98.00%
MW	36.62 kDa (predicted); 105.13 kDa (reducing conditions)
Uniprot ID	NP_004381.2
Expiration Date	6 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)