

Product Datasheet

Complement C5 Protein, Cynomolgus, Recombinant (His) (orb1953794)

Description	Complement C5 Protein, Cynomolgus, Recombinant (His) is expressed in HEK293 mammalian cells with His tag. The predicted molecular weight is 188.1 kDa and the accession number is XP_005580972.1.
Storage	-20°C
Tag	C-His
Note	For research use only
Application notes	A Certificate of Analysis (CoA) containing reconstitution instructions is included with the products. Please refer to the CoA for detailed information.
Purity	98.00%
MW	188.1 kDa (predicted); 106.3 and 75.8 kDa (reducing conditions)
Uniprot ID	XP_005580972.1
Expiration Date	6 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)