

Product Datasheet

PRTN3 Protein, Human, Recombinant (S203A, His) (orb1954101)

| | |
|--------------------------|---|
| Description | PRTN3 Protein, Human, Recombinant (S203A, His) is expressed in HEK293 mammalian cells with His tag. ... |
| Storage | -20°C |
| Tag | C-His |
| Note | For research use only |
| Application notes | A Certificate of Analysis (CoA) containing reconstitution instructions is included with the products. Please refer to the CoA for detailed information. |
| Purity | 98.00% |
| MW | 26.57 kDa (predicted) |
| Uniprot ID | P24158 |
| Expiration Date | 6 months from date of receipt. |