

## Product Datasheet

# Latent TGF beta 1 Protein, Human, Recombinant (His), Biotinylated (orb1954345)

<b>Description</b>	Latent TGF beta 1 Protein, Human, Recombinant (His), Biotinylated is expressed in HEK293 mammalian cells with His tag. The predicted molecular weight is 42.4 kDa and the accession number is A0A499FJK2.
<b>Storage</b>	-20°C
<b>Tag</b>	C-His
<b>Note</b>	For research use only
<b>Application notes</b>	A Certificate of Analysis (CoA) containing reconstitution instructions is included with the products. Please refer to the CoA for detailed information.
<b>Purity</b>	98.00%
<b>MW</b>	42.4 kDa (predicted)
<b>Uniprot ID</b>	<b>P01137</b>
<b>Expiration Date</b>	6 months from date of receipt.

---

**Biorbyt Ltd.**

7 Signet Court, Swann's Road,  
Cambridge, CB5 8LA, United Kingdom  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

**Biorbyt LLC.**

68 TW Alexander Drive,  
Durham, NC, 27713, United States  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)