

Product Datasheet

Noggin/NOG Protein, Human, Recombinant (orb1954625)

Description	Noggin/NOG Protein, Human, Recombinant is expressed in HEK293 mammalian cells. The predicted molecular weight is 23.1 kDa and the accession number is Q13253.
Storage	-20°C
Tag	Tag Free
Note	For research use only
Application notes	A Certificate of Analysis (CoA) containing reconstitution instructions is included with the products. Please refer to the CoA for detailed information.
Purity	SDS-PAGE: 99.1%; SEC-HPLC: 99.3%
MW	23.1 kDa (predicted); 30 kDa (reducing condition, due to glycosylation)
Uniprot ID	Q13253
Expiration Date	6 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)