

Product Datasheet

SCF Protein, Mouse, Recombinant (His), Biotinylated (orb1954669)

Description	SCF Protein, Mouse, Recombinant (His), Biotinylated is expressed in Baculovirus insect cells with His tag. The predicted molecular weight is 19.7 kDa and the accession number is P20826-1.
Storage	-20°C
Tag	C-His
Note	For research use only
Application notes	A Certificate of Analysis (CoA) containing reconstitution instructions is included with the products. Please refer to the CoA for detailed information.
Purity	98.00%
MW	19.7 kDa (predicted)
Uniprot ID	P20826
Expiration Date	6 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)