

## Product Datasheet

# RIOK1 Protein, Human, Recombinant (His & GST) (orb1954988)

<b>Description</b>	RIOK1, also known as RIO kinase 1, is a member of the RIO family of atypical serine protein kinases ...
<b>Storage</b>	-20°C
<b>Tag</b>	His, GST
<b>Note</b>	For research use only
<b>Application notes</b>	A Certificate of Analysis (CoA) containing reconstitution instructions is included with the products. Please refer to the CoA for detailed information.
<b>Purity</b>	98.00%
<b>MW</b>	93.4 kDa (predicted)
<b>Uniprot ID</b>	<a href="#">Q9BRS2</a>
<b>Expiration Date</b>	6 months from date of receipt.