

Product Datasheet

EphA4 Protein, Human, Recombinant (aa 570-986, His & GST) (orb1955260)

Description	EphA4 Protein, Human, Recombinant (aa 570-986, His & GST) is expressed in Baculovirus insect cells w...
Storage	-20°C
Tag	His, GST
Note	For research use only
Application notes	A Certificate of Analysis (CoA) containing reconstitution instructions is included with the products. Please refer to the CoA for detailed information.
Purity	98.00%
MW	75 kDa (predicted)
Uniprot ID	P54764
Expiration Date	6 months from date of receipt.