



Product Datasheet

CD11b antibody (orb19554)


Description	Goat polyclonal antibody to CD11b
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC
Target	ITGAM / CD11B
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:32000.WB: Approx 170kDa band observed in Human Brain (Cerebellum) lysates (calculated MW of 127kDa according to NP_000623.2). The observed molecular weight corresponds to earlier findings with different antibodies from other sources. Recommended concentration: 1-3 µg/ml.
Protein Sequence	RNDGEDSYRTQ
Clonality	Polyclonal
MW	127
Entrez	3684
Dilution Range	ELISA: 1:32000, IHC-P: 10 µg/ml
Expiration Date	12 months from date of receipt.



orb19554 Negative Control showing staini...



orb19554 (2µg/ml) staining of paraffin ...



Immunofluorescence analysis of U937 cell...