

## Product Datasheet

# Rad6/UBE2A Protein, Human, Recombinant (His) (orb1957186)

<b>Description</b>	Ubiquitin-conjugating enzyme E2 A (also known as HHR6A or UBE2A), encoded by human DNA repair genes ...
<b>Storage</b>	-20°C
<b>Tag</b>	N-His
<b>Note</b>	For research use only
<b>Application notes</b>	A Certificate of Analysis (CoA) containing reconstitution instructions is included with the products. Please refer to the CoA for detailed information.
<b>Purity</b>	98.00%
<b>MW</b>	19.2 kDa (predicted)
<b>Uniprot ID</b>	<b>P49459</b>
<b>Expiration Date</b>	6 months from date of receipt.