

## Product Datasheet

# EphA4 Protein, Human, Recombinant (His & hFc) (orb1957599)

<b>Description</b>	EphA4 Protein, Human, Recombinant (His & hFc) is expressed in HEK293 mammalian cells with His and hF...
<b>Storage</b>	-20°C
<b>Tag</b>	His, hFc
<b>Note</b>	For research use only
<b>Application notes</b>	A Certificate of Analysis (CoA) containing reconstitution instructions is included with the products. Please refer to the CoA for detailed information.
<b>Purity</b>	98.00%
<b>MW</b>	86.5 kDa (predicted)
<b>Uniprot ID</b>	<a href="#">A1MA61</a>
<b>Expiration Date</b>	6 months from date of receipt.