

Product Datasheet

EphA4 Protein, Human, Recombinant (His) (orb1958156)

Description	EphA4 Protein, Human, Recombinant (His) is expressed in HEK293 mammalian cells with His tag. The pre...
Storage	-20°C
Tag	C-His
Note	For research use only
Application notes	A Certificate of Analysis (CoA) containing reconstitution instructions is included with the products. Please refer to the CoA for detailed information.
Purity	98.00%
MW	60 kDa (predicted)
Uniprot ID	A1MA61
Expiration Date	6 months from date of receipt.