



Product Datasheet

PCBP4 antibody (orb19596)

Description	Goat polyclonal antibody to PCBP4
Species/Host	Goat
Reactivity	Canine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	PCBP4
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:128000.WB: Approx. 40kDa band observed in lysates of cell line HeLa (calculated MW of 41.4kDa according to NP_127502.1 and NP_127503.1). In transfected HEK293 transiently expressing PCBP4 a band of approx. 45kDa is observed. This band is not observed in the non-transfected HEK293. Recommended concentration: 1-3 µg/ml.
Protein Sequence	NGSKKAERQKFSPY
Clonality	Polyclonal
MW	43.9, 37.9, 41.5, 41.5
Entrez	57060
Dilution Range	ELISA: 1:128000, WB: 1-3 µg/ml
Expiration Date	12 months from date of receipt.



orb19596
(0.03ug/ml)
staining of
HeLa ce...



HEK293
overexpressing
PCBP4
(RC200749)
a...