

Product Datasheet

GFRAL/GFR alpha-like Protein, Mouse, Recombinant (His & Avi), Biotinylated (orb1959870)

Description	GFR alpha -like (GDNF receptor-alpha-like) is a distant member of the GDNFR family of proteins. Mature human GFR alpha-like is a 376 amino acid (aa) type I transmembrane protein. It contains a 333 aa extracellular domain, a 20 aa transmembrane domain and a 23 aa cytoplasmic domain. GFRAL is a brainstem-restricted receptor for GDF15 which regulates food intake, energy expenditure and body weight in response to metabolic and toxin-induced stresses. GFRAL/GFR alpha-like Protein, Mouse, Recombinant (His & Avi), Biotinylated is expressed in HEK293 mammalian cells with N-His-Avi tag. The predicted molecular weight is 40.1 kDa and the accession number is Q6SJE0-1.
Storage	-20°C
Tag	N-His-Avi
Note	For research use only
Application notes	Reconstitute the lyophilized protein in distilled water. The product concentration should not be less than 100 µg/ml. Before opening, centrifuge the tube to collect powder at the bottom. After adding the reconstitution buffer, avoid vortexing or pipetting for mixing.
Purity	98.00%
MW	40.1 kDa (predicted). Due to glycosylation, the protein migrates to 48-60 kDa based on Tris-Bis PAGE result.
Uniprot ID	Q6SJE0
Expiration Date	6 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)