



Product Datasheet

SLC34A2 antibody (orb195991)


Description	SLC34A2 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human SLC34A2. Please contact us for the exact immunogen sequence. The peptide is available as orb374834.
Target	SLC34A2
Form/Appearance	10 mM PBS, 0.02% sodium azide
Concentration	- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Polyclonal antibodies are purified by peptide affinity chromatography
MW	77 kDa
Uniprot ID	Q9DBP0 , Q9JJ09 , O95436
NCBI	001171470 , 001171469
Entrez	10568
Dilution Range	WB: 1:200, IHC-P: 1:300, IF/ICC: 1:300
Expiration Date	12 months from date of receipt.



WB analysis of Mouse small intestine (lane 1)



Western blot analysis of Rat liver (lane 1)



Mouse small intestine (lane 4), Rat kidn...