

Product Datasheet

DKK1 N terminal Domain Protein, Human, Recombinant (hFc & Avi), Biotinylated (orb1960721)

Description	Dickkopf-1 (Dkk1), the founding and best-studied member of the Dkk family, functions as an antagonist of canonical Wnt/ β -catenin. Dkk1 is considered to play a broad role in a variety of biological processes. DKK1 N terminal Domain Protein, Human, Recombinant (hFc & Avi), Biotinylated is expressed in HEK293 mammalian cells with C-hFc-Avi tag. The predicted molecular weight is 40.25 kDa and the accession number is O94907.
Storage	-20°C
Tag	C-hFc-Avi
Note	For research use only
Application notes	Reconstitute the lyophilized protein in distilled water. The product concentration should not be less than 100 μ g/ml. Before opening, centrifuge the tube to collect powder at the bottom. After adding the reconstitution buffer, avoid vortexing or pipetting for mixing.
Purity	98.00%
MW	40.25 kDa (predicted). Due to glycosylation, the protein migrates to 50-60 kDa based on Tris-Bis PAGE result.
Uniprot ID	O94907
Expiration Date	6 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)