

Product Datasheet

Neuropilin-1 Protein, Human, Recombinant (aa 22-644, Avi & His), Biotinylated (orb1961893)

| | |
|------------------------|---|
| Description | Neuropilin-1 (Npn-1, previously neuropilin; also CD304) is a 130 - 140 kDa type I transmembrane (TM)... |
| Storage | -20°C |
| Tag | C-Avi-6xHis |
| Note | For research use only |
| MW | 85-100 KDa (reducing condition) |
| Expiration Date | 6 months from date of receipt. |