

Product Datasheet

EFNA4 Protein, Human, Recombinant (hFc & His) (orb1962131)

Description	Ephrin-A4 is a member of the ephrin ligand family which binds members of the Eph receptor family. Al...
Storage	-20°C
Tag	C-hFc-6xHis
Note	For research use only
Application notes	Reconstitute the lyophilized protein in distilled water. The product concentration should not be less than 100 µg/ml. Before opening, centrifuge the tube to collect powder at the bottom. After adding the reconstitution buffer, avoid vortexing or pipetting for mixing.
Protein Sequence	Leu26-Gly171
Purity	98.00%
MW	45 KDa (reducing condition)
Uniprot ID	P52798
Expiration Date	6 months from date of receipt.