



Product Datasheet

Human IGF1R/Igf1 Receptor ELISA Kit (orb196294)

Description	Human IGF1R/Igf1 Receptor ELISA Kit
Reactivity	Human
Range	156 pg/ml - 10,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: E31-N932
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Application notes	Dilution ratio of 1:1, concentration in serum and plasma is around 500 pg/ml
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	154793 MW
Uniprot ID	P08069
Sensitivity	10 pg/ml
Expiration Date	Please enquire.



A typical standard curve for the Human I...



Human IGF1R PicoKine ELISA Kit standard ...