


Product Datasheet

Human IL-17E/IL-25 ELISA Kit (orb196298)

Description	Human IL-17E/IL-25 ELISA Kit
Reactivity	Human
Range	62.5 pg/ml - 4,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: Y33-G177
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Application notes	For serum, dilution ratio is 1:1 and concentration is less than the lowest standard, 62.5 pg/ml. For plasma, dilution ratio is 1:1 and concentration is 90 pg/ml
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	20330 MW
Uniprot ID	Q9H293
Sensitivity	10 pg/ml
Expiration Date	Please enquire.



Human IL-17E/IL-25 PicoKine ELISA Kit st...