


Product Datasheet

Human SPARCL1/SPARC-like protein 1 ELISA Kit (orb196299)

| | |
|----------------------------|---|
| Description | Human SPARCL1/SPARC-like protein 1 ELISA Kit |
| Reactivity | Human |
| Range | 46.9 pg/ml - 3,000 pg/ml |
| Tested Applications | ELISA |
| Immunogen | Expression system for standard: NS0; Immunogen sequence: I17-F664 |
| Storage | Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.) |
| Note | For research use only |
| Application notes | Dilution ratio of 1:200, concentration in serum and plasma is around 100 ng/ml |
| Sample Types | cell culture supernatants and serum. |
| MW | 75208 MW |
| Uniprot ID | Q14515 |
| Sensitivity | 10 pg/ml |
| Expiration Date | Please enquire. |



Human SPARCL1 PicoKine ELISA Kit standar...