


Product Datasheet

Human Mesothelin ELISA Kit (orb196302)

Description	Human Mesothelin ELISA Kit
Reactivity	Human
Range	156 pg/ml - 10,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: E296-G580
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Application notes	Dilution ratio of 1:10, concentration in serum and plasma is around 7.7-48.24 ng/ml
Sample Types	cell culture supernatants, serum, plasma (heparin, EDTA) , saliva and urine.
MW	68986 MW
Uniprot ID	Q13421
Sensitivity	10 pg/ml
Expiration Date	Please enquire.



Human Mesothelin PicoKine ELISA Kit stan...